

# Medical Vacuum Plant - British Standard



Precision UK's CPX Vacuum Plant is manufactured under BS EN 13485 Medical Devices: Quality Management Systems. It provides a centralised source of vacuum, available in simplex to multi-pump configurations dependant on individual requirement, to include complete package plant or modular format.

## Specific Benefits Include

- Fully complies to HTM 02-01, HTM 2022, C11, BS EN 7396, CE0086 Class IIb
- CE marked under Medical Device Directive 93/42/EEC
- CPX® branded solid British design
- Control and highly visible indication panel
- All interconnecting pipework is included
- Outlet pipework for connecting to the pipeline
- Anti-vibration mounts on all pumps
- Isolator to each pump unit as standard
- Alarm and BMS outputs for remote status indication
- Digital display available on request
- All settings are factory set and verified with certificates of inspection
- Each plant has a full warranty of 12 months from its commission date
- Options – Touchscreen, BMS Modbus, Digital Outputs, SMS Texting Alarms, Service Reminders



Packaged Plant