


Medical Air Plant - British Standard



Touchscreen Control Medical Air Dryer Unit

Precision UK's CPX Medical Air Plant is manufactured under BS EN 13485 Medical Devices: Quality Management Systems. It provides a centralised source of medical air, available in simplex to multi-pump configurations dependant on individual requirement.



Triplex Air Plant – HTM 02-01

Specific Benefits Include

- Fully complies to HTM 02-01, HTM 2022, C11, BS EN 7396, **CE0086 Class IIb**
- **CE marked under Medical Device Directive 93/42/EEC**
- CPX® branded solid British design
- Control and highly visible indication panel
- All interconnecting pipework is included
- Outlet pipework for connecting to the pipeline
- Anti-vibration mounts on all compressors
- Isolator to each pump unit as standard
- Alarm and BMS outputs for remote status indication
- **Digital display touchscreen with user interface**
- 400Kpa, 700Kpa and 1000Kpa output versions available
- All settings are factory set and verified with certificates of inspection
- Each plant has a full warranty of 12 months from its commission date